

EXPERIMENT No.: 1

TO PREPARE 100ml OF 0.1N SODIUM HYDROXIDE SOLUTION.

OBSERVATIONS:

No.	Initial Burette Reading	Final Burette Reading	DIFFERENCE (Volume used)
1.	0	7	7
2.	0	6.8	6.8
3.	0	6.2	6.2
Mean Volume Used =			6.77

CALCULATIONS:

$$N_1 V_1 = N_2 V_2$$

$$0.1 \times 10 = N_2 \times 100$$

$$\frac{0.1 \times 10}{100} = N_2$$

$$N_2 = 0.150$$

Now,

$$N_2 V_2 = N_1 V_1$$

$$0.150 \times V_2 = 0.1 \times 100$$

$$V_2 = \frac{0.1 \times 100}{0.150}$$

$$V_2 = 66.67$$

Now,

$$V_2 = 100 - 66.67$$

$$V_2 = 33.33 \text{ ml}$$

Required.

TO PREPARE 100ml OF 0.1N ACETIC ACID.OBSERVATIONS:

No.	Initial Burette Reading	Final Burette Reading	DIFFERENCE (Volume used)
1.	0	10.5	10.5
2.	10.5	20.5	10
3.	20	31.3	10.3
Mean Volume Used = 10.3 ml			

CALCULATIONS:

Normality of acid = $N_1 = ?$

Volume of Acid = $V_1 = 10 \text{ ml}$

Normality of Base = $N_2 = 0.1 \text{ N}$

Volume of Base = $V_2 = 10.3 \text{ ml}$

$$N_1 V_1 = N_2 V_2$$

$$N_1 \times 10 \text{ ml} = 0.1 \text{ N} \times 10.3$$

$$N_1 = \frac{0.1 \text{ N} \times 10.3}{10}$$

$N_2 = 0.103 \text{ ml}$ (but required normality is 0.1)

So, $N_2 = N_3 = \text{stock sol. prepared}$

Normality required = $N_R = 0.1 \text{ N}$

Volume required = $V_R = 100 \text{ ml}$

Normality of stock = $N_S = 0.103 \text{ ml}$

Vol. of stock = $V_S = ?$

So,

$$N_S V_S = N_R V_R$$

$$(0.103)(V_S) = (0.1)(100)$$

$$V_S = \frac{(0.1)(100)}{0.103}$$

$$V_S = 97 \text{ ml}$$

→ To make 100ml of 0.1N
add water to 97ml

So, water added

$$= 100 - V_S$$

$$= 100 - 97 = 3 \text{ ml water}$$

EXPERIMENT No.: _____

Date: _____

TO PREPARE 100ml OF 0.1N HYDROCHLORIC ACID.OBSERVATIONS:

No.	Initial Burette Reading	Final Burette Reading	DIFFERENCE (Volume used)
1.	0	15	15
2.	15	30	15
3.	30	45	15
Mean Volume Used =			15 ml

CALCULATIONS:Normality of Stock HCl = $N_1 = ?$ Volume of Stock HCl = $V_1 = 10 \text{ ml}$ Normality of standard NaOH = $N_2 = 0.1 \text{ N}$ Volume of Standard NaOH = $V_2 = 15 \text{ ml}$ So, $N_1 V_1 = N_2 V_2$

$$N_1 \times 10 = 0.1 \times 15$$

$$N = \frac{1.5}{10} = 0.15 \text{ N}$$

So Normality of Stock HCl is 0.15 N which is more than required normality i.e. 0.1 N

Normality of req. sol = $N_R = 0.1 \text{ N}$ Volume " " " = $V_R = 100 \text{ ml}$ Normality of Stock sol = $N_S = 0.15$ Volume " " " = $V_S = ?$

8

$$N_S V_S = N_R V_R$$

$$0.15 \times V_S = 0.1 \times 100$$

$$V_S = \frac{10}{0.15}$$

$$= 66.66 \text{ ml}$$

To make 100ml sol
add H₂O to
stock solWater to be
added = $100 - V_S$
= $100 - 66.66$
= 33.34 ml


DETECTION OF GLUCOSE IN A GIVEN SOLUTION.

APPARATUS: Test tubes, Test tube holder, Test tube Rack, Pipettes, Beaker, Water bath, Burner, Glass slide, cover slip and Microscope.

REAGENTS: Molisch's Reagent, Iodine Reagent, Benedict's Reagent, Barfoed's Reagent, Seliwanoff's Reagent, Osazone Mixture, Conc. HCl, Conc. H₂SO₄, 5% NaOH & Original solution (Glucose solution).

No.	TESTS	OBSERVATION	INFERENCE
1	MOLISCH'S TEST	Reddish violet ring	Carbohydrate
		at the junction of two liquids	confirmed
2	IODINE TEST	No colour change	Polysaccharide
		occur	absent
3	BENEDICT'S TEST	Brick red colour	Presence of red-
			sugar confirmed
4	BARFOED'S TEST	formation of red-ppt	Monosaccharide
		before 7 min	confirmed
5	SELIWANOFF'S TEST	No colour change	Aldo group
		noticed	confirmed
6	OSAZONE TEST	Yellow, needle shaped	Glucosazone
		crystals	confirmed

SKETCH OF THE OSAZONE CRYSTALS



RESULT: Glucose confirmed

Date: _____

EXPERIMENT No. _____


DETECTION OF FRUCTOSE IN A GIVEN SOLUTION.

APPARATUS: Test tubes, Test tube holder, Test tube Rack, Pipettes, Beaker, Water bath, Burner, Glass slide, cover slip and Microscope.

REAGENTS: Molisch's Reagent, Iodine Reagent, Benedict's Reagent, Barfoed's Reagent, Seliwanoff's Reagent, Osazone Mixture, Conc. HCl, Conc. H₂SO₄, 5% NaOH & Original solution (Fructose solution).

No.	TESTS	OBSERVATION	INFERENCE
1	MOLISCH'S TEST	Reddish violet ring appears at the junction of two liquids	Carbohydrate confirmed
2	IODINE TEST	No colour change	Absence of polysaccharide or disaccharide
3	BENEDICT'S TEST	Red brick colour formed on heating	Reducing sugar present
4	BARFOED'S TEST	Formation of red ppt before 7 min	Monosaccharide present
5	SELIWANOFF'S TEST	Formation of red colour precipitate observed	Keto group present
6	OSAZONE TEST	Yellow needle shaped crystals observed	Fructose present

SKETCH OF THE OSAZONE CRYSTALS




RESULT: Fructose present

DETECTION OF GALACTOSE IN A GIVEN SOLUTION.

APPARATUS: Test tubes, Test tube holder, Test tube Rack, Pipettes, Beaker, Water bath, Burner, Glass slide, cover slip and Microscope.

REAGENTS: Molisch's Reagent, Iodine Reagent, Benedict's Reagent, Barfoed's Reagent, Seliwanoff's Reagent, Osazone Mixture, Conc. HCl, Conc. H₂SO₄, 5% NaOH & Original solution (Galactose solution).

No.	TESTS	OBSERVATION	INFERENCE
1	MOLISCH'S TEST	Reddish violet ring	Carbohydrate
		observed	present
2	IODINE TEST	No colour change	Polysaccharide
		observed	absent
3	BENEDICT'S TEST	Red-brick colour	Reducing sugar
		formed on heating	present
4	BARFOED'S TEST	Formation of red ppt	Monosaccharide
		before 7 min	present
5	SELIWANOFF'S TEST	No red ppt / colour	Aldo group
		formed	present
6	OSAZONE TEST	Fluffy ball shaped	Galactosazone
		crystals formed	Formed

SKETCH OF THE OSAZONE CRYSTALS

RESULT: Galactose present

EXPERIMENT No. _____


DETECTION OF LACTOSE IN A GIVEN SOLUTION.

APPARATUS: Test tubes, Test tube holder, Test tube Rack, Pipettes, Beaker, Water bath, Burner, Glass slide, cover slip and Microscope.

REAGENTS: Molisch's Reagent, Iodine Reagent, Benedict's Reagent, Barfoed's Reagent, Seliwanoff's Reagent, Osazone Mixture, Conc. HCl, Conc. H₂SO₄, 5% NaOH & Original solution (Lactose solution).

No.	TESTS	OBSERVATION	INFERENCE
1	MOLISCH'S TEST	Reddish violet ring formed	Carbohydrate present
		No blue colouration	Polysaccharide absent
2	IODINE TEST	Green coloured sol	Reducing sugar present
		formed on heating	Disaccharide present
3	BENEDICT'S TEST	Formation of red ppt after 7 min	Aldo group present
		Red colour not formed	present
4	BARFOED'S TEST	Red colour not formed	Aldo group present
		formed	present
5	SELIWANOFF'S TEST	Red colour not formed	Aldo group present
		formed	present
6	OSAZONE TEST	Fluffy ball shaped crystals formed	Lactosazone present
		crystals formed	present

SKETCH OF THE OSAZONE CRYSTALS



RESULT: Lactose present

EXPERIMENT No. _____

Date: _____

DETECTION OF MALTOSE IN A GIVEN SOLUTION.


APPARATUS: Test tubes, Test tube holder, Test tube Rack, Pipettes, Beaker, Water bath, Burner, Glass slide, cover slip and Microscope.

REAGENTS: Molisch's Reagent, Iodine Reagent, Benedict's Reagent, Barfoed's Reagent, Seliwanoff's Reagent, Osazone Mixture, Conc. HCl, Conc. H₂SO₄, 5% NaOH & Original solution (Maltose solution).

No.	TESTS	OBSERVATION	INFERENCE
1	MOLISCH'S TEST	Red violet ring	Carbohydrate
		formed	present
2	IODINE TEST	No blue colouration	Polysaccharide
		observed	absent
3	BENEDICT'S TEST	Olive-green colour	Reducing sugar
		observed on heating	present
4	BARFOED'S TEST	Formation of red	Disaccharide
		ppt after 7 min	present-
5	SELIWANOFF'S TEST	Red colour not	Aldo group
		formed	present
6	OSAZONE TEST	Sunflower shaped	Maltose
		crystals formed	present-

SKETCH OF THE OSAZONE CRYSTALS

RESULT: Maltose present



EXPERIMENT No. _____


DETECTION OF SUCROSE IN A GIVEN SOLUTION.

APPARATUS: Test tubes, Test tube holder, Test tube Rack, Pipettes, Beaker, Water bath, Burner, Glass slide, cover slip and Microscope.

REAGENTS: Molisch's Reagent, Iodine Reagent, Benedict's Reagent, Barfoed's Reagent, Seliwanoff's Reagent, Osazone Mixture, Conc. HCl, Conc. H₂SO₄, 5% NaOH & Original solution (Sucrose solution).

No.	TESTS	OBSERVATION	INFERENCE
1	MOLISCH'S TEST	Formation of red	Carbohydrate present
		violet ring	
2	IODINE TEST	No colour change	Polysaccharide absent
		observed	
3	BENEDICT'S TEST	No colour change	Non-reducing sugar present
		observed	
4	BARFOED'S TEST	-	-
5	SELIWANOFF'S TEST	Cherry red coloured	keto group present
		ppt's formed	
6	PERFORM THE FOLLOWING TESTS AFTER HYDROLYSIS.		
a	BENEDICT'S TEST	Brick red colour	Red. sugar present
		formed	
b	BARFOED'S TEST	Red ppt. observed	Monosaccharide present
		before 7. mins	
c	OSAZONE TEST	Needle shaped	C ₇ glucose / Fructose present
		crystals are observed	

SKETCH OF THE OSAZONE CRYSTALS



RESULT: Sucrose present

EXPERIMENT No. _____


DETECTION OF STARCH IN A GIVEN SOLUTION.

APPARATUS: Test tubes, Test tube holder, Test tube Rack, Pipettes, Beaker, Water bath, Burner, Glass slide, cover slip and Microscope.

REAGENTS: Molisch's Reagent, Iodine Reagent, Benedict's Reagent, Barfoed's Reagent, Seliwanoff's Reagent, Osazone Mixture, Conc. HCl, Conc. H₂SO₄, 5% NaOH & Original solution (Starch solution).

No.	TESTS	OBSERVATION	INFERENCE
1	MOLISCH'S TEST	Reddish violet ring	Carbohydrate present
		is formed	
2	IODINE TEST	Colour turns	Polysaccharide present
		blue	
3	BENEDICT'S TEST	No change in colour	Red sugar absent
		of solution	
4	BARFOED'S TEST	No ppt formed	Polysaccharide present
		upon prolonged heating	
5	SELIWANOFF'S TEST	No change in	keto sugar absent
		colour	
6	PERFORM THE FOLLOWING TESTS AFTER HYDROLYSIS.		
a	BENEDICT'S TEST	Brick red ppt	Red-sugar present
		formed	
b	BARFOED'S TEST	Red ppt formed	Monosaccharide present
		before 7 min	

SKETCH OF THE OSAZONE CRYSTALS



RESULT: Starch confirmed

Date: _____

EXPERIMENT No. _____


DETECTION OF UNKNOWN CARBOHYDRATE IN A GIVEN SOLUTION.

APPARATUS: Test tubes, Test tube holder, Test tube Rack, Pipettes, Beaker, Water bath, Burner, Glass slide, cover slip and Microscope.

REAGENTS: Molisch's Reagent, Iodine Reagent, Benedict's Reagent, Barfoed's Reagent, Seliwanoff's Reagent, Osazone Mixture, Conc. HCl, Conc. H₂SO₄, 5% NaOH & Original solution (Unknown Carbohydrate solution).

No.	TESTS	OBSERVATION	INFERENCE
1	MOLISCH'S TEST	Reddish violet ring formed	Carb present.
2	IODINE TEST	No colour change	Polysaccharide absent.
3	BENEDICT'S TEST	Red brick ppt formed	Red sugar present.
4	BARFOED'S TEST	Red ppt formed after 7 min	Disaccharide present.
5	SELIWANOFF'S TEST	No red colour formed	Keto group absent.
6	PERFORM THE FOLLOWING TESTS AFTER HYDROLYSIS OF THE UNKNOWN SOLUTION (IF REQUIRED)		
A	BENEDICT'S TEST	-	-
B	BARFOED'S TEST	-	-
C	OSAZONE TEST	Fluffy ball shaped crystals formed	Lactose present.

SKETCH OF THE OSAZONE CRYSTALS



RESULT: Lactose present